




Piston Nebulizer

Model: NE100



- 2in1 cup for easy cleaning, assembly and disassembly
- MMAD $\leq 2.6 \mu\text{m}$
Fine Particle Dose (FPD): 80-85%
- Neb Rate $\geq 0.35 \text{ ml/min}$
- Consistent fine particle size for efficient respiratory treatment
- Powerful piston compressor
- Mouthpiece and masks for adult and child included
- Economic design



Efficient Treatment



Hygienic Design
Easy to Clean



Adult and Child
Masks



Patent Nos :
ZL 2013 3 0528746.6
ZL 2015 3 0231205.6
Other patents pending

Nebulizer Kit
Consisting of:
• 2in1 cup (model N4)
• Mouthpiece



Aerosol performance tested
with 2.5% NaF solution by Cascade Impactor

Cascade Impactor


MMAD	FPD
2.591 μm	81.69%

FPD: Fine particle dose, the percentage of particle size less than 5.0 μm

Hygienic "2in1 Cup" Design



Particle Size Distribution



Model	Qty per carton	Carton volume
NE100	6 pcs	0.039 cbm/ctn
QTY per 20'	QTY per 40'	QTY per 40'HQ
4308 pcs	8616 pcs	9996 pcs

● Storage and Transportation Condition
-20°C to 60°C (-4°F to 140°F)
10-90% RH, 700-1060hPa

www.rossmax.com
OBM_Epaper_NE100_N4_ver1702