




Power & Efficiency Suction Unit

Model: V5

Triple Protection consist of

- 1 Overflow valve 2 Antibacterial filter 3 Fool-proof design



Patent Nos:
ZL 2016 3 0298084.1
Other patents pending



It's Powerful

- High performance pump delivers -600mmHg vacuum and 30 liter/min flow rate

It's Protected "Triple Protection"

- Overflow valve plus antibacterial filter and fool-proof design that prevents back flow of fluid and airborne contaminants from entering the pump

It's Convenient

- Easy to store using tubing holder, tubing depository and accessory compartment

It's Ergonomic

- Comprehensive dashboard with ergonomic view angle ensures a smooth operation
- Easy to open with button operated lid open design

Technical Specification

- Suction class: High Vacuum / High Flow
 - Max vacuum @ pump inlet: -80 kPa/ -600 mmHg
 - Max air flow @ pump inlet: 30 L/min
 - Power source: AC100-120V, 50/60Hz or AC220-240V, 50/60Hz
 - Collection canister capacity: 1400ml (Polycarbonate, PC)
 - Power Consumption: ≤ 55VA
- * Subject to technical modification without prior notice

Model	Qty per carton	Carton volume
V5	2 pcs	0.070 cbm/ctn
QTY per 20'	QTY per 40'	QTY per 40'HQ
480 pcs	960 pcs	1120 pcs
QTY per pallet	48 pcs	